

# THERAPOD SPECIALIST SEATING



## Before

When seated for long periods on a standard chair, the pelvis rotates and the spinal column bends forward, resulting in a tendency to slouch. The epressure placed on the spine by this posture is an inevitable precursor of back pain.



## After

All Status Seating solutions promote a natural forward tilt, allowing the spine to maintain a healthy 'S' curve. System Therapod® chairs also uniquely enable the user to customise lumbar support to their specific body shape.

## THERAPOD 24/7

- Extra high back Therapod shell with detachable neck support
- Therapod patented backcare system with adjustable strapping to shape the chair back for comfort and support
- V rated flexible polyurethane back foam
- Upholstered in 24/7 fabric as standard
- Steel reinforced armrest bolted to both seat back and seat ensures massive structural strength
- High density pressure sensitive seat foam
- Upgraded 30 stone gas lift
- 28 inch heavy duty polished aluminium base with heavy duty castors
- Heavy duty two lever rocker action fitted with extra strong tension spring



### Therapod 24/7 measurements

Seat depth to back	480
Seat width maximum	500
Seat maximum height	560
Back width at maximum	500
Overall height of chair from ground maximum	1420
Gap between arms	510

4500 £650.00

## THERAPOD 7470

The Therapod 7470 enables the larger person to enjoy smart, modern office seating but with the comfort of a chair designed to accommodate their physical needs. Like all Therapods this model employs System Therapod® adjustable strapping to enable the user to shape the seat back. The model 7470 has an internal seat width of 600mm, a full 100mm bigger than the standard chair width, which enables well-built employees to fit easily between the seat arms without embarrassment.

To accommodate the additional weight the chair features:

- A heavy duty 28 inch polished aluminium base with 50mm hooded chrome castors
- An upgraded gas lift to accommodate up to 190kg (thirty stone)
- Extra thick (100mm) pressure sensitive seat foam
- Heavy duty British manufactured advanced synchronised mechanism with up-rated tension spring
- Extra heavy duty T arm with quick release height adjustment
- Guaranteed for two years under normal useage



### Therapod 7470 measurements

Seat depth to back	480
Seat width maximum	600
Seat maximum height	550
Back width at maximum	500
Overall height of chair from ground maximum	1290
Gap between arms	580


7470 £589.00

Inside the back of every Therapod chair is a flexible polycarbonate/ABS shell mounted on to a steel ribcage for added strength. The shell is made to flex and contort to fit the characteristics of the users back by the adjustment of four hidden support straps. Known as dynamic seating, this unique system allows lateral movement and twisting actions while providing optimal spinal support.

Status uses pressure sensitive seat foam on all its products. This unique product responds to body heat and shapes itself to spread the weight evenly so as to equalise the pressure over the whole surface. This provides even greater comfort and stability for the user.



## CUSTOMISING YOUR THERAPOD



1. To customise the Therapod to your back, first unzip the back cover to expose the adjustable straps.
2. Loosen off the four adjustable straps by unfastening the velcro tabs.
3. Pull out each strap to the required tension and then velcro it securely to retain the setting. Support increases by the tension of each strap.
4. Re-zip the chair cover and enjoy the long-term comfort. Once a chair is set up for its user, it senses and responds to that individual's needs whatever seating position they adopt.

### THERAPOD 5250

- High back Therapod shell
- Quick lift height adjustment giving 84mm of travel over seven positions
- Therapod patented backcare system with adjustable strapping to shape the chair back for comfort and support
- V rated flexible polyurethane back foam
- Multi functional arms, height adjustable with armpads that swivel and move backwards and forwards
- Pressure sensitive seat foam
- 150kg (23 stone) Level S rated gas lift
- 27 inch polished aluminium base
- 50mm twin wheel castors
- British manufactured advanced synchronised action featuring:
  - Patented fingertip control side tension adjustment with eight incremental settings in a 270° turn
  - Anti kick back safety feature as required by EN1335 (1999)
  - Seven incremental back rest angle lock settings
  - 23° back rest recline
  - 16° of seat movement starting at 4° forward
  - Minimal shirt and knee lift for superior comfort



### Therapod 5250 measurements

Seat depth to back	460
Seat width maximum	510
Seat maximum height	550
Back width at maximum	500
Overall height of chair from ground maximum	1240
Gap between arms	470

**5250** £499.00

### THERAPOD 5500

- Extra high back Therapod shell with detachable neck support
- Quick lift height adjustment giving 84mm of travel over seven positions
- Therapod patented backcare system with adjustable strapping to shape the chair back for comfort and support.
- V rated flexible polyurethane back foam
- Multi-functional arms, height adjustable with armpads that swivel and move backwards and forwards.
- Pressure sensitive seat foam
- 150kg (23 stone) Level S rated gas lift
- 27 inch polished aluminium base
- 50mm twin wheel castors
- British manufactured advanced synchronised action featuring:
  - Patented fingertip control side tension adjustment with eight incremental setting in a 270° turn
  - Anti kick back safety feature as required by EN1335 (1999)
  - Seven incremental back rest angle lock settings
  - 23° back rest incline
  - 16° of seat movement starting at 4° forward
  - Minimal shirt and knee lift for superior comfort



### Therapod 5500 measurements

Seat depth to back	460
Seat width maximum	500
Seat maximum height	550
Back width at maximum	500
Overall height of chair from ground maximum	1410
Gap between arms	470

**5500** £545.00